PROBENTEILER

SAMPLE DIVIDER

PROBENTEILER

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Sample Divider

SICON Sample Dividing Equipment is used to handle bulk materials samples of any description in the laboratories or industry, research and development establishments, etc., anywhere the world over. This equipment has been supplied in a great variety of versions to meet individual requirements and /or for service in complete sampling systems.

SICON Sample Dividing Equipment impresses by its heavy-duty design and resultant sampling accuracy, the number one prerequisite for meaningful analytical results.

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ROTARY SAMPLE DIVIDER MODEL PT

The Rotary Sample Divider divides larger sample quantities and extracts samples from flowing bulk material.

Depending upon model size these dividers will sample and divide all flowable materials up to a particle size of 60 mm, continuous automatic sample division ensures that the fraction sampled will be highly representative of the whole. The unit essentially comprises the following parts:

Motor, v-belt drive, stainless steel dividing tube, top cone, inside cone, sample outlet pipe, reject outlet and guard.

The dividing tube is rotated by a motor and V-belt. The material is charged through a bore in the drive pulley and flows through the revolving dividing tube into the sample outlet and against the undercone. The reject material passes through the undercone bore. A lever permits to open and adjust the sample outlet. Depending upon the width of this opening, the following division ratio can be defined:

$$\text{Sampling ratio} = \frac{\text{Division circle } x}{\text{Sample slot width } y}$$

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
Rotary sample dividers (Model K) are supplied as complete units with Vibratory Feeders and funnel.

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